**GIMMI**®

**Surgical Instruments**

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D - 91054 Erlangen  
roman.carbon@chir.imed.uni-erlangen.de

## Abdomen and Thorax Ring System acc. Carbon

Due to the different age and weight of the patient in pediatric surgery there are limitations to the universal usage of passive holding systems such as retractors and rings. Especially in the premature and neonatal surgery there remains a need for adaptable holding units as the patients are so tiny and the OP-team has only a small area to work on. The "hook holding assistant", if one is needed, cannot usually see into the operation field, and due to this his hook can be somewhat ineffective at times. Because the structure of these tiny patients are so delicate, dangerous hook movements can be the result.

A new modulating ring system, specially for premature and neonatal pediatric thorax and abdomen surgery, has now been developed. This ring system accommodates the delicate constitution of the rib cage wall and torso and the extreme structure elasticity of babies who only weighs about 500 - 2500g.

Contrary to the presently used diminutive ring and valve systems, this lightly oval shaped Carbon neonatal ring system is made of stainless steel of only 1 mm thickness and has an inside diameter of 10 x 8,5 cm or 13 x 10 cm (see illustration 1 and 2). This polished ring is fully flashless and has soft contours. Both the ring and valves are made of stainless instrument steel.

illustration 1

illustration 2

The small ring only weighs about 20 g and as one valve only weighs 7 g, the total weight of six valves is 65 g. For better fixation the valves carry the ring on an outer edge with a wavy atraumatic grind whereby every square offers four fixation variations.

The valves are similar to atraumatic lid hooks or similar to Langenbeck retractors (see illustration 2) and they themselves offer three fixation position possibilities for the width of the OP-field and can be held at the fine shaped end and manipulated accordingly. (illustration 2, 3, 4)


illustration 3

## Abdomen and Thorax Ring System acc. Carbon

The smallest OP field is a diagonal of 2,5 cm (illustration 3) and the largest has a dimension of 6 cm (illustration 4). Due to the wavy form of the ring, an adequate number of valves can be used simultaneously. The use of four valves usually results in a diamond or trapez shaped OP field, an additional two valves then gives an hexagonal shaped field and if eight valves are used the OP field becomes nearly circular in shape. These different modelling variations are important as, due to the high elasticity of the integuments of a premature baby, the tissue between the valves can become taut and this can block the view to the OP field. By using the Langenbeck type of valves there is the possibility to partly push aside the parenchyma organs such as liver and lungs or the hollow organs such as stomach intestines or bladder.


illustration 4

Six valves seem to be the ideal configuration, e.g. for medial laparotomy (see illustration 5, 6) whereby initially, after the ring has been put into place, two valves are placed opposite each other, thus already providing extensive exposure.

Finally, two valves can then be attached cranially and caudally at a 120° angle to each other and thus protect the cranial structures (navel vein, lig. falciform) and caudal structures (urachus, bladder). Eversion of the intestine is harmless as the ring and valves are fully flashless and polished and, due to their "soft" form, no squashing is possible. During lateral thoraxotomy, after the ring has been put into place, the distension of the tissue can proceed: cranial and caudal ribs are separated with two valves. Further valves can be used to pull the tissue and soft parts away from the rib cage.

illustration 5

illustration 6

### Case example:

Premature baby, 1400 g, strangulation ileus through ductus omphaloentericus. Use of Carbon ring system in pediatric surgery. Six valves are used to form a hexagonal shaped OP field for better exposure and all valves are fixed into the first position. Eversion of the intestine and the laying of organs onto the holding system is harmless due to the fact that the surface of the holding system is polished and has a "soft" shape.

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

Valve, **small**, for abdominal and thoracic retractor acc. Carbon

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Valve, **large**, for abdominal and thoracic retractor acc. Carbon

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Valve, **extra large**, for abdominal and thoracic retractor acc. Carbon

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Valve, **round**, for abdominal and thoracic retractor acc. Carbon

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Dissecting Scissors, curved  
with microtome cut, 10,0 cm

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Vessel Forceps, curved, 9,0 cm  
very soft tension

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Neonatal-Forceps, curved, 15,0 cm  
very soft tension

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We recommend the following SET:

### Basic Instruments, modified acc. Carbon:

<b>D.9299.90</b>	Frame for abdominal and thoracic retractor acc. CARBON	- 1 Stück
<b>D.9299.99</b>	Frame for abdominal and thoracic retractor acc. CARBON	- 1 Stück
<b>D.9299.96</b>	Valve, small, for abdominal and thoracic retractor acc. CARBON	- 2 Stück
<b>D.9299.97</b>	Valve, large, for abdominal and thoracic retractor acc. CARBON	- 2 Stück
<b>D.9299.98</b>	Valve, round, for abdominal and thoracic retractor acc. CARBON	- 6 Stück
<b>AB.1189.18</b>	GERALD Dissecting Forceps, 17,5 cm,	- 3 Stück
<b>B.9934.18</b>	GERALD Tissue Forceps, 1x2 teeth, 17,5 cm	- 3 Stück
<b>C.0015.090S</b>	HARTMANN Vessel Forceps, curved, 9,0 cm	- 6 Stück
<b>A.9011.10</b>	Dissecting Scissors, curved, 10,0 cm, with mikrotom cut	- 3 Stück
<b>I.5595.01B</b>	Neonatal-Forceps acc. CARBON, curved, 15,0 cm	- 2 Stück
<b>R.9030.20</b>	Container, 296 x 275 x 115 mm	- 1 Stück (o. Abb.)
<b>R.9030.21</b>	Wrapping drapes, Cotton	- 1 Stück (o. Abb.)
<b>R.9030.22</b>	Sterilizing tray, 60 mm	- 1 Stück (o. Abb.)
<b>R.9030.23</b>	Silicon mats	- 1 Stück (o. Abb.)
<b>R.9030.24</b>	Identification label	- 1 Stück (o. Abb.)

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

**GIMMI® GmbH**

Surgical Instruments  
Endoscopic Technology

Carl-Zeiss-Straße 6  
D-78532 Tuttlingen

Telefon ++ 49 74 61 / 96 590-0

Telefax ++ 49 74 61 / 96 590-33

email [info@gimmi.de](mailto:info@gimmi.de)

website [www.gimmi.de](http://www.gimmi.de)

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Ltd. Oberarzt an der Abteilung für Kinderchirurgie